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Technology & Platform

- Sample Management
 - Samportal
 - DNA extraction
- SNP Genotyping
 - Fine mapping
 Sequenom MassARRAY System
 - Human whole-genome scan
 Affymetrix GeneChip
 Affymetrix GeneTitan
- STRP Genotyping
- Next Generation Sequencing
 - illumina HiSeq 2000
 - illumina MiSeq
- BioIT Services: iQuestion, GenePipe, HiPipe Basic, HiPipe Professional, Varpipe
- Statistical Analysis: GWAS

Human Whole-Genome Genotyping using Affymetrix GeneChip System / Affymetrix GeneTitan System

Affymetrix GeneChip System

Affymetrix GeneTitan System

Affymetrix Human Whole-Genome Genotyping Arrays (Statistical analysis included)

	ASI	CHB	TWB	SNP 6.0
Technology platform	Affymetrix GeneTitan	Affymetrix GeneTitan	Affymetrix GeneTitan	Affymetrix GeneChip
Population	Asian	Han Chinese in Beijing, China	Han Chinese	all
SNP probes	600,307	642,832	653,291	906,600
DNA required	300 ng (15 ng/uL*20 uL)	300 ng (15 ng/uL*20 uL)	300 ng (15 ng/uL*20 uL)	500 ng (50 ng/uL*10 uL)
Average call rate	99.63%	99.54%	99.69%	98.80%
Method	Ligation	Ligation	Ligation	Hybridization
CNV probes	none	none	none	946,000
Throughput (samples/week)	285	285	285	36
DNA source	cell lines, blood, and saliva	cell lines, blood, and saliva	cell lines, blood, and saliva	cell lines and blood
Service code	C2-15	C2-16	C2-42	C2-14

APPLICATIONS

SNP Genotyping using Human-Whole Genome Genotyping Arrays

- Linkage analysis
- Whole-genome association study
- Population genetics
- Chromosomal copy number analysis
- Cancer genetics
- Linkage disequilibrium

Candidate gene/loci SNP Genotyping using Sequenom MassARRAY System

- **MALDI-TOF (Matrix-Assisted Laser Desorption / Ionization Time-Of-Flight) Mass Spectrometry based platform**
- **Multiplex assay design: cost effective solution**

APPLICATIONS

SNP Genotyping using Sequenom MassARRAY System


- High-throughput genotyping of customer selected SNPs.
- Fine mapping of SNPs in candidate genes or pre-localized regions for both human and mouse samples.
- Whole-genome SNP genotyping for disease gene mapping in ENU study and speed congenics mouse model establishment.

ABI 3730 DNA Analyzer

APPLICATIONS

STRP Genotyping

Journal of Food, Agriculture & Environment (2014) 12 (2):541-553

- Whole genome linkage analysis
- PCR fragment analysis
- Kinship analysis: most commonly in paternity testing
- Forensic identification
- Genetic improvement (strain improvement)
- Animal breeding
- Examination of genetic diversity and phylogeny
- Conservation of biodiversity

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	illumina HiSeq 2000	illumina MiSeq				
Sequencing Chemistry	Sequencing by synthesis (Reversible terminator)					
Amplification	Bridge amplification					
Read Length	100 bp x 2 (Paired-end)	1x50 bp	2x150 bp*	2x250 bp**	2x300 bp	150 bp
No. of Reads/ Flow Cell	175-200 million / lane (8 lanes/Flow Cell)	12~15 million			20~25 million	
Output/run	480-600 Gb / 2 Flow cells	600~750 Mb / FC	3.5~4.5 Gb / FC	6~7.5 Gb / FC	12~15 Gb / FC	3~3.8 Gb / FC
Multiplex	12 barcodes	~ 96 barcodes				
Instrument Run Time/run	13 days	1 day	1 day	2 days	3 days	1 day
Data Processing /run	2.5 - 4 days (depends on project)	2.5 - 4 days (depends on project)				

* 2x150 bp Nano: 0.8~1 million (240~300 Mb / FC); Micro: 3~4 million (1~1.2 Gb)

** 2x250 bp Nano: 0.8~1 million (400~500 Mb / FC)

APPLICATIONS
Next Generation Sequencing

illumina HiSeq 2000

- Genome re-sequencing
- *de novo* sequencing
 for small genome
- Transcriptome sequencing
 - mRNA, stranded
 - rRNA-depleted, stranded
- Exome sequencing
- ChIP sequencing
- RNA IP sequencing
- Ribosome profiling
- Small RNA sequencing
- Methylation sequencing (coming soon)

illumina MiSeq

- *de novo* sequencing for small genome
- Whole genome re-sequencing for small genome
- Metagenomics (16s rRNA sequencing)
- Targeted re-sequencing (two-step PCR amplification, gene panels, Nextera XT tagmentation)
- Small RNA sequencing

Next Generation RNA Sequencing

Application	Transcriptome	Small RNA	RNA IP
Sample source	Total RNA	Total RNA	RNA IP
Sample requirement-Amount and Conc.	① Agilent Bioanalyzer Amount: 100ng ~ 2μg Conc.>50 ng/μL RNA Integrity Number (RIN) > 8	① Agilent Bioanalyzer Amount: >3μg Conc.>100 ng/μL RNA Integrity Number (RIN) > 8	① Agilent Bioanalyzer Amount:>10ng Size> 100 nt
Sample requirement-OD ratio	② NanoDrop OD260/280: 1.8~2.0 OD260/230: 1.8~2.2		
Sample requirement-others	③ DNase treated ④ No contaminants such as organic solvent, inhibitors, etc. ⑤ In RNase-Free Water (DEPC Water)		
Service Code for library prep	1) C2-36 mRNA sequencing sample preparation 2) C2-37 rRNA-Depleted RNA sequencing sample preparation	C2-45 Small RNA sequencing sample preparation	C2-36 mRNA sequencing sample preparation
Platform	HiSeq	HiSeq or MiSeq	HiSeq or (MiSeq)

Next Generation DNA Sequencing

Application	Genome Resequencing	Exome Sequencing	CHIP sequencing	Targeted Resequencing	
Sample source	Genomic DNA	Genomic DNA	CHIP DNA	Amplicon	Amplicon
Sample requirement-Amount and Conc.	① Qubit Amount: 100ng~1μg Conc.>100 ng/μL	① Qubit Amount: > 3 μg Conc.>100 ng/μL	① Qubit & Agilent Bioanalyzer Amount: > 10 ng Size: 150~500 bp	① Qubit Amount: 100ng~1μg Conc.>100 ng/μL	① Qubit Amount: > 10 ng Conc.>15 ng/μL
Sample requirement-OD ratio	② NanoDrop OD260/280: 1.8~2.0 OD260/230: >1.7				
Sample requirement-others	③ RNase treated ④ No degradation of gDNA (examined by 1% agarose gel electrophoresis) ⑤ No contaminants such as organic solvent, inhibitors, etc. ⑥ In molecular grade water or 10 mM Tris-HCl (pH 7.5–8.5)				
Service Code for library prep	1) C2-27 Human whole genome sequencing, 30X 2) C2-29 DNA sequencing sample preparation	C2-48 SureSelect Exome Enrichment (Human/Mouse)	C2-50 ChIP sequencing sample preparation	C2-29 DNA sequencing sample preparation	C2-51 DNA library sample preparation by Nextra XT
Platform	HiSeq or MiSeq	HiSeq	HiSeq	MiSeq	MiSeq

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Transcriptome	Small RNA	RNA IP	Genome Resequencing	Exome Sequencing	Targeted Resequencing	CHIP sequencing
Differential Expression Analysis Gene fusion detection variant detection De novo assembly	Index trimming, QC, Grouping (read-count format), Alignment, Filter (rRNA, tRNA, mRNA, sRNA, genome), Expression Level, isomiR identification, Novel miRNA identification, Differential Expression	peak calling peak annotation	DNA basic analysis: mapping, Variant calling, annotation of SNVs and indels (Gene information, location in coding, UTR, or intron regions, splice site change, nonconservative/ conservative amino acid substitution, related known disease, affected protein domain, predicted splicing regulation), novel SNVs/indels detection / annotation, and exome CNV detection. Advanced analysis: customized grouping analysis. LOH analysis Amplicon and whole genome CNV detection Low frequency mutation detection			peak calling motif search peak annotation

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

Human whole genome genotyping


GWAS analysis

Validation Fine mapping Replication


Joint analysis


SCIENCE nature genetics The NEW ENGLAND JOURNAL of MEDICINE

Causal variant analysis


Variants and genes interpretation

GenePipe V3.0

NGS data analysis


HiPipe

High-throughput candidate gene sequencing

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Stage I
Sample Collection
- Electronic data capture system
- Sample management system

Stage II
Screening & Discovery
- Genotyping technology
- Genomics technology
- Bioinformatics service
- Biostatistics service

Stage III
Validation
- Genotyping technology
- Genomics technology

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
Total solutions for genotyping, NGS and NGS/GWAS data analysis!

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